User Manual

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	试剂盒组成	规格	二抗	
一抗	Mouse IgG1, κ Isotype Control Antibody (同型对照)	20 μL	FITC/PE Goat Anti-mouse IgG	
	Purified Anti-rat CD90	20 μL		
	Purified Anti-rat CD44	20 μL		
	Purified Anti-rat CD73	20 μL		
	Purified Anti-rat CD34	20 μL		
	Purified Anti-rat CD45	20 μL		
	Mouse IgG2a, κ Isotype Control Antibody (同型对照)	20 μL		
	Purified Anti-rat CD11b/c	20 μL		
	Armenian Hamster IgG Isotype Control Antibody (同型对照)	20 μL		FITC/PE Goat Anti-hamster IgG
	Purified Anti-mouse/rat CD29	20 μL		
二抗	FITC Goat Anti-mouse IgG Antibody	80 μL		
	FITC Goat Anti-hamster IgG Antibody	20 μL		
	PE Goat Anti-mouse IgG Antibody	80 μL		
	PE Goat Anti-hamster IgG Antibody	20 μL		

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1. Prockop DJ. Marrow stromal cells as stem cells for nonhematopoietic issues. *Science* 1997;276(5309):71_4.
2. Pittenger MF, Mackay AM, Beck SC, Jaiswal RK, Douglas R, osca JD, et al. Multilineage potential of adult human mesenchymal stem cells. *Science* 1999;284(5411):143_7.
3. Krampera M, Pizzolo G, Aprili G, Franchini M. Mesenchymal stem cells for bone, cartilage, tendon and skeletal muscle repair. *Bone* 2006;39(4):678_83.
4. Ohnishi S, Nagaya N. Prepare cells to repair the heart: mesenchymal stem cells for the treatment of heart failure. *Am Nephrol* 2007;27(3):301_7.
5. Dezawa M, Hoshino M, Nabeshima Y, Ide C. Marrow stromal cells: implications in health and disease in the nervous system. *Curr Mol Med* 2005;5(7):723_32.
6. Humphreys BD, Bonventre JV. Mesenchymal stem cells in acute kidney injury. *Ann Rev Med* 008;59(1):311_25.
7. Giordano A, Galderisi U, Marino IR. From the laboratory bench to the patients beside: an update on clinical trials with mesenchymal stem cells. *J Cell Physiol* 2007;211(1):27_35.
8. Lee OK, Coathup MJ, Goodship AE, Blunn GW. Use of mesenchymal stem cells to facilitate bone regeneration in normal and chemotherapy-treated rats. *Tissue Eng* 2005;11(11_12):1727_35.
9. Peng L, Jia Z, Yin X, et al. Comparative analysis of mesenchymal stem cells from bone marrow cartilage, and adipose tissue. *StemCells Devel* 2008;17(4):761_74.

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