



## User Manual

# ***OriCell™ Strain C57BL/6 Mouse Mesenchymal Stem Cells with RFP (MSCs/RFP)***

**Cat. No. MUBMX-01201**

# Table of Contents

<b>Contents and Storage</b> .....	<b>3</b>
<b>Product Introduction</b> .....	<b>3</b>
<b>Cell Characteristics and Identity</b> .....	<b>3</b>
<b>Product Application</b> .....	<b>4</b>
<b>General Handling Principles</b> .....	<b>4</b>
<b>Culturing OriCell™ Strain C57BL/6 Mouse MSCs/RFP</b>	
<b>Thawing and Establishing OriCell™ Strain C57BL/6 Mouse MSCs/RFP</b> .....	<b>4</b>
<b>Passaging Cyagen OriCell™ Strain C57BL/6 Mouse MSCs/RFP</b> .....	<b>6</b>
<b>Differentiation of OriCell™ Strain C57BL/6 Mouse MSCs/RFP</b> .....	<b>7</b>
<b>Cryopreservation of OriCell™ Strain C57BL/6 Mouse MSCs/RFP</b> .....	<b>11</b>
<b>Appendix</b> .....	<b>13</b>
<b>Troubleshooting</b> .....	<b>13</b>
<b>Related Products</b> .....	<b>14</b>
<b>References</b> .....	<b>14</b>
<b>Technical Support</b> .....	<b>15</b>

## CONTENTS AND STORAGE

Product Name	Strain C57BL/6 Mouse Mesenchymal Stem Cells with RFP
Catalog No.	MUBMX-01201
Amount per Vial	1×10 <sup>6</sup> Cells
Cryopreserved At	Ninth Passage
Storage Condition	Liquid Nitrogen



**CAUTION:** Please handle this product as a potentially biohazardous material. This product contains dimethyl sulfoxide (DMSO), a hazardous material, in the freezing medium.

## PRODUCT INTRODUCTION

Mesenchymal stem cells (MSCs) are multipotent stem cells that can differentiate into a variety of cell types including osteocytes, adipocytes, and chondrocytes. MSCs proliferate quickly and are capable of generating a local immunosuppressive microenvironment, thus contributing to their wide application potentials in tissue engineering, cell therapy, and gene therapy.

OriCell™ C57BL/6 Mesenchymal Stem Cells with RFP are derived from the bone marrow of C57BL/6 mice, cultured as monolayer, and then have been transfected with a lentiviral construct containing a RFP expression motif. They have a strong capacity for self-renewal while maintaining their multipotency.

In addition, these cells have been tested for:

- **Exogenous Factors:** bacterial/fungal contamination, mycoplasma contamination, and endotoxin contamination.
- **Characteristics:** post-thaw viability, cell cycle, verification of undifferentiated state, and differentiation potential.

This product is intended for laboratory research use only. It is not intended for diagnostic, therapeutic, clinical, household, or any other applications.

## CELL CHARACTERISTICS AND IDENTITY

- Strong capacity to expand. Can be passaged at least 5 times.
- Multipotent differentiation ability along the osteogenic, chondrogenic, and adipogenic lineages.
- Positive for CD29, CD44 and Sca-1 (> 70%), and negative for CD117 and CD31 (< 5%) in flow cytometry assays.

## PRODUCT APPLICATIONS

Strain C57BL/6 Mouse MSCs/RFP have become a popular research target due to their potential use in regenerative medicine and tissue engineering (in areas such as cardiovascular, neural, and orthopedic disease).

OriCell™ Strain C57BL/6 Mouse MSCs/RFP can be used as cell models to evaluate the immunoreactions, proliferation, immigration, and differentiation of MSCs both *in vivo* and *in vitro*.

## GENERAL HANDLING PRINCIPLES

1. Aseptic handling of the product is necessary throughout.
2. Once the cells have been established, always freeze several vials of OriCell™ Strain C57BL/6 Mouse MSCs/RFP as a backup.



**Note:** The OriCell™ Strain C57BL/6 Mouse MSCs/RFP can be frozen/thawed at least two times.

3. For general maintenance of cells, we recommend the seeding density to be 2.0-3.0 × 10<sup>4</sup> cells/cm<sup>2</sup>.
4. For general maintenance of cells, we recommend that the medium is changed if it becomes acidic (the pH indicator in the medium appears yellow). In general, change the growth medium every three days.
5. Do not let OriCell™ C57BL/6 MSCs/RFP overgrow as it will result in contact inhibition. When the cells are 80-90% confluent, subculturing the cells is strongly recommended.



**Note:** We strongly recommend the use of OriCell™ culture media and other related reagents for optimal results.

## THAWING AND ESTABLISHING OriCell™ Strain C57BL/6 Mouse MSCs/RFP

### Materials Required

- OriCell™ Mouse Mesenchymal Stem Cell Growth Medium (Cat. No. MUXMX-90011)

### Thawing and Establishing Strain C57BL/6 Mouse MSCs/RFP

1. Pre-warm the fully supplemented (complete) OriCell™ Mouse MSC Growth Medium to 37°C.
2. Add 9 mL of OriCell™ Mouse MSC Growth Medium to a 15 mL conical tube.
3. Remove the cryovial of OriCell™ Strain C57BL/6 Mouse MSCs/RFP from liquid nitrogen.
4. Quickly thaw the cryovial in a 37°C water bath until the last ice crystal disappears. For optimal results, be sure to finish the thawing procedure within 3 minutes. Be careful not to submerge the entire vial. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.



**Note:** Results will be less than optimal if the cells are thawed for more than 3 minutes.

5. As soon as the cells are completely thawed, disinfect the outside of the cryovial with 70% v/v ethanol.
6. Use a pipette to transfer the cells to the 15 mL conical tube containing OriCell™ Mouse MSC Growth Medium inside a biosafety cabinet. Be careful not to introduce any bubbles during the transfer process.
7. Rinse the vial with 1 mL of the medium to reduce cell loss. Subsequently transfer this 1 mL of cell suspension into the conical tube.
8. Gently mix the cell suspension by slowly pipetting up and down. Be careful not to introduce any bubbles.
9. Centrifuge the cell suspension at 250 x g for 5 minutes.
10. Carefully aspirate off as much of the supernatant as possible and add 2-3 mL of fresh OriCell™ Mouse MSC Growth Medium (pre-warmed to 37°C).
11. Gently resuspend the cells in OriCell™ Mouse MSC Growth Medium.
12. Seed the cells into a T25 flask and add a sufficient amount of OriCell™ Mouse MSC Growth Medium. Gently rock the culture flask to evenly distribute the cells.
13. Incubate the flask at 37°C inside a 5% CO<sub>2</sub> humidified incubator.
14. The next day, change the medium with fresh growth medium (pre-warmed to 37°C).
15. Change the growth medium every three days thereafter.
16. When the cells are approximately 80-90% confluent, they can be dissociated with Trypsin-EDTA and passaged.



**Note: Changing Medium**

1. Warm an appropriate amount of medium to 37°C in a sterile container. Replace the spent medium with the pre-warmed, fresh medium. Once completed, return the flask to the incubator.
2. Avoid repeated warming and cooling of the medium. If the entire content is not needed for a single procedure, transfer only the required volume to a sterile secondary container.



**Fig. 1.** OriCell™ Strain C57BL/6 Mouse Mesenchymal Stem Cells with RFP are established.

## **PASSAGING OriCell™ Strain C57BL/6 Mouse MSCs/RFP**

### **Materials Required**

- 0.25% Trypsin-0.04% EDTA (Cat. No. TEDTA-10001)
- Phosphate-Buffered Saline (1×PBS) (Cat. No. PBS-10001)
- OriCell™ Strain C57BL/6 Mouse Mesenchymal Stem Cells with RFP (Cat. No. MUBMX-01201)
- OriCell™ Mouse Mesenchymal Stem Cell Growth Medium (Cat. No. MUXMX-90011)

### **Passaging OriCell™ Strain C57BL/6 Mouse MSCs/RFP**

1. Pre-warm the OriCell™ Mouse MSC Growth Medium, 1×PBS, and 0.25% Trypsin-0.04% EDTA solution to 37°C.
2. Carefully aspirate the spent medium from the 80-90% confluent monolayer of MSCs.
3. Add 1×PBS (6 mL for T75 flask, 3 mL for T25 flask). Be careful not to disturb the monolayer. Gently rock the flask back and forth to rinse the monolayer.
4. Aspirate 1×PBS off and discard.
5. Repeat steps 3-4 two or three times.
6. Add 0.25% Trypsin-0.04% EDTA solution (2-3 mL for T75 flask, 1 mL for T25 flask). Gently rock the flask back and forth to ensure that the entire monolayer is covered with the Trypsin-EDTA solution. Allow trypsinization to continue until the majority of the cells (approximately 80%) are rounded up. At this point, gently tap the side of the flask to release the majority of cells from the culture flask surface.



**Important:** Avoid leaving cells exposed to the trypsin longer than necessary (no more than two minutes if using Cyagen's trypsin-EDTA solution). Care should also be taken that the cells are not forced to detach prematurely as this may result in clumping.

7. After the cells are visibly detached, immediately add the pre-warmed OriCell™ Mouse MSC Growth Medium (6 mL for T75 flask, 3 mL for T25 flask) to neutralize the trypsinization.
8. Gently pipette the medium over the cells to dislodge and resuspend the cells. Repeat 5-6 times until all the cells are dissociated from the flask and evenly dispersed into a single cell suspension.
9. Transfer the dissociated cells into a 15 mL conical tube.
10. Centrifuge at 250 x *g* for 5 minutes.
11. Carefully aspirate off as much of the supernatant as possible.
12. Add 2 mL of OriCell™ Mouse MSC Growth Medium to the conical tube and gently resuspend the cells thoroughly.
13. Plate the cells into appropriate flasks. OriCell™ Strain C57BL/6 Mouse MSCs/RFP can be split at 1:2 or other appropriate ratios.
14. Add an appropriate amount of medium to the cells. Incubate the cells at 37°C inside a 5% CO<sub>2</sub> humidified incubator.



**Note:** Care should be taken to avoid introducing bubbles during pipetting.

## Additional Tips

### Time to Change Medium

It is recommended to change the culture medium if there are too many dead cells after passaging.

It is recommended to change the culture medium whenever the medium becomes acidic, even if the cells do not reach 80-90% confluency. The pH indicator in the culture medium will appear yellow when acidic.

### Time to Subculture

When OriCell™ Strain C57BL/6 Mouse MSCs/RFP are 80-90% confluent, it is recommended that the cells be subcultured. Do not let the cells overgrow as it will result in contact inhibition.

## OriCell™ Strain C57BL/6 Mouse MSCs DIFFERENTIATION USING OriCell™ DIFFERENTIATION MEDIA

OriCell™ Strain C57BL/6 Mouse MSCs/RFP can differentiate into a variety of cell types including osteocytes, adipocytes, and chondrocytes.

### Osteogenic Differentiation

## Materials Required

OriCell™ Mesenchymal Stem Cell Osteogenic Differentiation Medium (Cat. No. GUXMX-90021)

## Osteogenesis Protocol



**Note:** The protocol listed below is for 6-well tissue culture plates.

1. Culture the OriCell™ Strain C57BL/6 Mouse MSCs/RFP in OriCell™ Mesenchymal Stem Cell Growth Medium at 37°C in a 5% CO<sub>2</sub> humidified incubator.
2. When cells are approximately 80-90% confluent, they can be dissociated with 0.25% Trypsin-0.04% EDTA (Cat. No. TEDTA-10001).
3. Re-seed the MSCs in the growth medium at  $3 \times 10^4$  cells/cm<sup>2</sup> in a 6-well tissue culture plate pre-coated with 0.1% gelatin solution.
4. Incubate the cells at 37°C inside a 5% CO<sub>2</sub> humidified incubator.
5. When cells are approximately 60-70% confluent, carefully aspirate off the growth medium from each well and add 2 mL of OriCell™ Mesenchymal Stem Cell Osteogenic Differentiation Medium.
6. Feed cells every three days for 2-4 weeks by completely replacing the medium with fresh OriCell™ Mesenchymal Stem Cell Osteogenic Differentiation Medium (pre-warmed to 37°C).
7. After 2-4 weeks of differentiation, cells can be fixed and stained with alizarin red S.




**Note:** To prevent osteoblasts from detaching, it is recommended to change half of the medium every two days before analysis.

## Alizarin Red S Staining Analysis

1. After the cells have differentiated, remove the osteogenic differentiation medium from the wells and rinse with 1x phosphate-buffered saline (PBS). Fix cells with 2 mL of 4% formaldehyde solution for 30 minutes.
2. Rinse wells twice with 1x PBS. Stain the cells with 1 mL alizarin red S working solution for 3-5 minutes.
3. Rinse wells 2-3 times with 1x PBS.
4. Cells can now be visualized and analyzed under a microscope.





**Fig. 3** OriCell™ Strain C57BL/6 Mouse MSCs/RFP are differentiated into osteocytes and are stained with alizarin red S.

## Adipogenic Differentiation

### Materials Required

OriCell™ Mesenchymal Stem Cell Adipogenic Differentiation Medium (Cat. No. GUXMX-90031)

### Adipogenesis Protocol




**Note:** The protocol listed below is for 6-well tissue culture plates.

1. Culture the OriCell™ Strain C57BL/6 Mouse MSCs/RFP in the OriCell™ Mesenchymal Stem Cell Growth Medium at 37°C in a 5% CO<sub>2</sub> humidified incubator.
2. When cells are approximately 80-90% confluent, they can be dissociated with 0.25% Trypsin-0.04% (Cat. No. TEDTA-1000).
3. Reseed the MSCs in growth medium at 2x10<sup>4</sup> cells/cm<sup>2</sup> in a 6-well tissue culture plate with a medium volume of 2 mL per well.
4. Incubate the cells at 37°C in a 5% CO<sub>2</sub> humidified incubator.
5. Feed the cells every three days until they are 100% confluent or post-confluent. Induction of adipogenic differentiation at post-confluency is strongly recommended.
6. When the cells are 100% confluent or post-confluent, carefully aspirate off the spent growth medium from the wells and add 2 mL of OriCell™ Mesenchymal Stem Cell Adipogenic Differentiation medium A (induction medium) per well.
7. Three days later, change the medium to OriCell™ Mesenchymal Stem Cell Adipogenic Differentiation medium B (maintenance medium) by completely replacing the spent medium A.
8. 24 hours later, change the medium back to MSC Adipogenic Differentiation medium A.
9. To optimally differentiate MSCs into adipogenic cells, repeat the cycle of induction and maintenance at least three times.
10. After three to five cycles of induction and maintenance, culture the cells in OriCell™ Mesenchymal Stem Cell Adipogenic Differentiation medium B for an

additional 4-7 days until the lipid droplets are big, round enough. During these days period, change the medium every three days.

### Oil Red O Stain Analysis

1. After the cells have differentiated, remove the MSC maintenance medium from the wells and rinse with 1x phosphate-buffered saline (PBS). Fix cells with 2 mL of 4% formaldehyde solution for 30 minutes.
2. Rinse wells twice with 1x PBS and stain cells with 1 mL of oil red O working solution (3:2 dilution with distilled water and filter with filter paper) for 30 minutes.
3. Rinse wells 2-3 times with 1x PBS.
4. Cells can now be visualized and analyzed under a microscope.



**Fig.4** OriCell™ Strain C57BL/6 Mouse MSCs/RFP are differentiated into adipocytes and are stained with oil red O.

## Chondrogenic Differentiation

### Materials Required

OriCell™ Mesenchymal Stem Cell Chondrogenic Differentiation Medium  
(Cat. No. GUXMX-90041)


### Chondrogenesis Protocol

1. Calculate the total number of MSC pellet cultures required for your experiment ( $2.5 \times 10^5$  MSCs are needed to form each chondrogenic pellet). Transfer this amount of cells into an appropriate culture tube.
2. Wash the MSCs with Incomplete Chondrogenic Medium. Centrifuge the cells at  $150 \times g$  for 5 minutes at room temperature and then aspirate off the supernatant. Resuspend the cells in 1 mL of Incomplete Chondrogenic Medium per  $7.5 \times 10^5$  cells. Centrifuge again at  $150 \times g$  for 5 minutes and then aspirate off the medium.
3. Resuspend the MSCs in Complete Chondrogenic medium to a concentration of  $5.0 \times 10^5$  cells/mL.

4. Aliquot 0.5 mL ( $2.5 \times 10^5$  cells) of the cell suspension into 15 mL polypropylene culture tubes. Centrifuge the cells at  $150 \times g$  for 5 minutes at room temperature. DO NOT aspirate the supernatant or resuspend the pellet.
5. Loosen the caps of the tubes one half turn in order to allow gas exchange, and incubate the tubes at  $37^\circ\text{C}$  in a humidified atmosphere of 5%  $\text{CO}_2$ . Do not disturb the pellets for 24 hours.
6. Feed the cell pellets every 2-3 days by completely replacing the medium in each tube (to avoid aspirating the pellets when aspirating the medium, attach a sterile 1-200  $\mu\text{L}$  pipette tip to the end of the aspirating pipette). Add 0.5 mL of freshly prepared Complete Chondrogenic Medium to each tube.
7. After replacing the medium, flick the bottom of the tube to ensure that the pellet is freefloating. Loosen the caps and return the tubes to the  $37^\circ\text{C}$  incubator.
8. Chondrogenic pellets should be harvested after 14-28 days in culture. Pellets may be formalin-fixed and paraffin-embedded for alcian blue stain analysis.

### Alcian Blue Staining Procedure

1. The tissue sample should be formalin-fixed and paraffin-embedded already.
2. Staining procedure:
  - a) Deparaffinize slides and hydrate to distilled water.
  - b) Stain in alcian blue solution for 30 minutes.
  - c) Wash in running tap water for 2 minutes.
  - d) Rinse in distilled water.
  - e) Visualize under a light microscope and capture images for analysis. Blue staining indicates synthesis of proteoglycans by chondrocytes.



**Fig.5** OriCell™ Strain C57BL/6 Mouse MSCs/RFP are differentiated into cartilages and are stained with alcian blue.

## CRYOPRESERVATION OF CELLS USING OriCell™CRYOPRESERVATION MEDIA

OriCell™NCR Protein-Free Cryopreservation Medium (Cat. No. NCPF-10001) is a protein-free, ready-to-use freezing medium. Itschemically-definedandprotein-freeformulationhas beenoptimizedtostemcellsandprimarycells,thusgreatlyenhancingtheviabilityand integrityofthesecellsbyprotectingthemfromdamageduringtheone-stepfreeze-thaw procedure.Unlikeotherconventionalfreezingmedia,whichrequireaslowprogrammedfreeze, thisproductallowsthecellstobedirectly frozen at -80°C.

### Cryopreservation



**Note:** Change the culture medium with fresh growth medium 24 hours before freezing.

1. Collectcellsthat are in the logarithmicgrowthphase.Perform a cell counttodeterminetheviablecell density.
2. Centrifugethecellsfor3-5minutesat250 x *g*and20°C.Removeanddiscardthe supernatantusingapipette.
3. Resuspendthecell pelletintheOriCell™NCRProtein-FreeCryopreservationMediumatacell density of10<sup>5</sup>-10<sup>6</sup>cells/mL.
4. Dispensealiquotsofthecell suspensionintocryogenicstoragevialsthatareproperly labeled.
5. Placethevialsdirectlyina -80°Cfreezer.After24hours,transferthefrozenvialstoliquid nitrogenforlong-term preservation.

## APPENDIX

### Troubleshooting

The table below lists some potential problems and solutions for culturing MSCs.

Problem	Cause	Solution
Low cell recovery rate	The storage condition does not meet the requirements	Purchase a replacement and store in liquid nitrogen for long-term preservation.
	Thawing of the cells takes too long	Thaw cells for no more than 3 minutes.
	Cells are incompletely recovered after thawing	After aspirating off medium, wash the tube with culture medium twice and transfer all of the cells to the dish.
	Cells are handled roughly	Care should be taken to avoid introducing bubbles during pipetting. Also avoid vortexing and high-speed centrifugation.
	Medium is not pre-warmed	Warm medium to 37°C before recovery.
Slow cell growth	Mycoplasma contamination	Discard the cells in question and disinfect the laboratory environment before recovering the next batch of cells.
	Over digestion	Wash the cells with PBS 2-3 times to remove serum prior to trypsinization (serum will inhibit the function of trypsin).
		Control the digestion time.
Plating density is too low	Increase the plating density.	
Cell aging	Inappropriate serum and medium	Use Cyagen tailor-made culture media. If other serum and media products are used, please perform validation to ensure compatibility.
	Dead cells are not removed promptly	Change the medium next day after recovery to ensure removal of all dead cells.
	Cell Contamination	Discard the cells in question and disinfect the laboratory environment before recovering the next batch of cells.
	Plating density is too low	Some stem cells can secrete factors to support cell growth. Therefore, a certain degree of plating density must be maintained; otherwise, it will lead to cell proliferation slow down and cell aging.

Cell aging	Over digestion	Wash the cells with PBS 2-3 times to remove serum prior to trypsinization (serum will inhibit the function of trypsin).
		Control the digestion time.
	The passaging time is not appropriate	The cells should be subcultured when reaching 80-90% confluency in order to avoid contact inhibition.
Cells show spontaneous differentiation	DMSO is not completely removed during cell recovery	Wash the cells with pre-warmed medium 2-3 times during recovery.
Ineffective induction of cell differentiation	Differentiation reagents need to be optimized	Use Cyagentailor-made differentiation media.
	Cell passage is too high	Use cells at a low original passage number.

## RELATED PRODUCTS

Product	Catalog Number
<b>OriCell™ Mouse Mesenchymal Stem Cell Growth Medium</b>	MUXMX-90011
<b>OriCell™ Mesenchymal Stem Cell Osteogenic Differentiation Medium</b>	GUXMX-90021
<b>OriCell™ Mesenchymal Stem Cell Adipogenic Differentiation Medium</b>	GUXMX-90031
<b>OriCell™ Mesenchymal Stem Cell Chondrogenic Differentiation Medium</b>	GUXMX-90041
<b>0.25% Trypsin-0.04% EDTA</b>	TEDTA-10001
<b>Phosphate-Buffered Saline (1xPBS)</b>	PBS-10001
<b>OriCell™ NCR Protein-Free Cryopreservation Medium</b>	NCPF-10001

## REFERENCES

Alexandra Peister, Jason A. Mellad, and Benjamin L. Larson. (2004) Adult stem cells from bone marrow (MSCs) isolated from different strains of inbred mice vary in surface epitopes, rates of proliferation, and differentiation potential. *BLOOD* 1: 1662-1668.

Masoud Soleimani, and Samad Nadn. (2009) A protocol for isolation and culture of mesenchymal stem cells from mouse bone marrow. *NETURE* 4: 102-106.

Philippe Tropel, Daniele Noel, and Nadine Platet. (2004) Isolation and characterisation of mesenchymal stem cells from adult mouse bone marrow. *Experimental Cell Research* 295: 395-406.

Lindolfo da Silva Meirelles, Pedro Cesar Chagastelles, and Nance Beyer Nardi. (2006) Mesenchymal stem cells reside in virtually all post-natal organs and tissues. *Journal of Cell Science* 119: 2204-2213.

**Cyagen Biosciences reserves all rights on the technical documents of its OriCell™ cell culture products. No part of this document may be reproduced or adapted for other purposes without written permission from Cyagen Biosciences.**