

## OriCell Mouse Embryonic Stem Cell Cryopreservation Medium

Catalog No. MUXES-07031

### Product Description:

OriCell Mouse Embryonic Stem Cell Cryopreservation Medium is a serum-containing cell freezing medium of Mouse Embryonic Stem Cells.

The product is intended for laboratory research use only, and not for clinical and house hold applications, or other uses.

### Features

- Special for Mouse Embryonic Stem Cells
- Prepare in Basal Medium of cells (DMEM)
- Formulated with Cell-Qualified Fetal Bovine Serum
- Formulated with DMSO(10%)

### Packing

50 mL/ Bottle

20 mL/ Bottle

### Instructions:

#### Cryopreservation

1. Select the cells in logarithmic growth phase according to the methods used to collect cells in vitro. Perform a cell count to determine the viable cell density.
2. Centrifuge the cells for 3-5minutes at 250g, 20°C. Remove the supernatant using a pipette.
3. Resuspend the precipitate with Cryopreservation Medium at the cell density of  $10^5$ - $10^6$  cells/mL.
4. Dispense aliquots of the cell suspension into cryogenic storage vials that are properly labeled.
5. Place the vials in a programmed controlled rate freezing apparatus and make sure the cooling rate of 1°C/min.
6. Transfer the frozen vials to liquid nitrogen for preservation after 24 hours.

## Thawing

1. Remove the cryovial of frozen cells from storage and quickly thaw the vial in 37°C water bath.
2. Gently dilute and mix the cells with 8mL-10mL of cell culture medium.
3. Centrifuge the cells for 3-5minutes at 250g, 20°C.
4. Carefully aspirate as much of the supernatant as possible and add 2-3mL of cell culture medium to resuspend the cells.
5. Transfer the cells into appropriate growth vessel and incubate at 37°C in a 5% CO<sub>2</sub> humidified incubator.

## Stability/Storage:

Store at 2-8°C. The product is stable at 2-8°C for up to 1 year. The product should be discarded beyond the labeled expiration date.

For optimal performance, repeated warming/cooling and freeze-thawing should be avoided.

## Quality Control:

Mouse Embryonic Stem Cell Cryopreservation Medium is performance tested on Mouse Embryonic Stem Cell.

Standard evaluation includes:

1. Sterility test (bacteria, fungi, mold and mycoplasma)
2. pH test
3. Osmolality
4. Endotoxin

**Cyagen Biosciences reserves all rights on the technical documents of its OriCell cell culture products. No part of this document may be reproduced or adapted for other purposes without written permission from Cyagen Biosciences.**